

## How to purchase a commercially fabricated splint.

1. An ErinoakKids' clinician will schedule assessment/measurements for a Benik or Professional Orthopedic Products (POPs) splint.
2. The clinician will complete and send the order form to the vendor on your behalf.
3. Provide your email address and confirm that we can share your email address, street address or phone number with our vendor.
4. The clinician will request a follow-up appointment for your fitting (~6-8 weeks after your initial appointment)
5. You will receive your quote via email **or** regular mail from the vendor. (depending on your preference).
6. Once you receive your quote you will need to make payment directly to the company i.e. Benik or POPs. (Please keep in mind quotes are only valid for 30 days)
  - Benik uses UPS for their delivery – this is a more expensive option, however it is a faster option (2-3 days transit time) and more predictable.
  - You may also choose to go with the US postal service with completion of delivery done by Canada Post. This will take longer in transit and is more **unpredictable**.
7. You will receive the splint at your home address in ~ 4-6 weeks.
8. Be sure to attend your follow up appointment to ensure that the splint fits properly. The splint may need to be modified (for example; thermoplastic modeled, neoprene trimmed, etc.) Please do not wear your splint until those modifications are in place.

### Contact Information:



#### Benik Corporation

11871 Silverdale Way NW #107  
Silverdale, WA 98383 USA

Call toll free: 1-800-442-8910

Phone: (360) 692-5601

Fax: (360) 692-5600

Email: [info@benik.com](mailto:info@benik.com)



#### Professional Orthopedic Products

5735 Ferrier  
Mont-Royal, Quebec, H4P 1N3

Toll free: 1-866-334-4441

Phone: 514 484-4441

Fax: 514-484-4431

Email: [csr@poportho.com](mailto:csr@poportho.com)